

HDROZONE™

Ultra-Low Temperature Sterilization (37–42 °C):

Polymer Safety, Lumen Performance, and Clinical Outcomes

Hasan Tahsin ÖZBEK • Bülent DEVECİ • Yüksel ERGÜN

White Paper - **05.11.2025**

ABSTRACT

In modern surgical and endoscopic practice, sterilization is not merely a process of microbial elimination; it is a critical process in terms of preserving device integrity and maintaining clinical safety on a sustainable basis.

Flexible endoscopes, robotic surgical instruments, and polymer-based medical devices present specific requirements during sterilization due to their sensitivity to heat and oxidative stress.

In this study, the HDROZONE™ approach has been evaluated on the basis of ultra-low temperature (37–42 °C) operation and controlled reactive chemistry principles, with polymer safety, intraluminal sterilization performance, and impact on clinical outcomes addressed systematically.



KEYS

HDROZONE, Polymer, Lumen, TASS, Residuals, Toxicity, Microstructure, Biofilm, Penetration, Stability

INTRODUCTION

Sterilization technologies are the primary determinant of the safe use of medical devices.

Conventional approaches are based on:


- High temperatures
- Aggressive chemical agents

However, medical devices in current use — flexible endoscopes, robotic surgical instruments, and polymer-based complex systems — are not fully compatible with classical sterilization methods.

This situation creates a critical need to balance sterilization efficacy with material integrity.

DEFINITION OF THE CLINICAL PROBLEM





The Endoscope and Polymer Era

Modern devices are no longer predominantly metal-based. They incorporate:

- PEEK
- PTFE
- Silicone
- Elastomers
- Optical fiber systems

These materials are:

- Sensitive to high temperatures
- Susceptible to microstructural deformation
- Vulnerable to chemical stress

Clinical Risks in the Field

Key problems observed in clinical practice:

- Uncertainty regarding intraluminal sterilization
- Endoscope damage (micro-cracks)
- Degradation of image quality
- Risk of Toxic Anterior Segment Syndrome (TASS)
- Reuse safety concerns

These conditions directly translate into clinical outcomes.

ULTRA-LOW TEMPERATURE (37–42 °C): THE CRITICAL THRESHOLD

Temperature Range of Current Systems

Many low-temperature systems operate in the 45–55 °C range. This range is close to the critical threshold values for polymers.

Thermal Behavior of Polymers

The key concept for polymers is the Glass Transition Temperature (T_g). As the T_g is approached:

- Chain mobility increases
- Microstructure becomes relaxed
- Deformation begins

Effect of the 37–42 °C Range

The HDROZONE™ operating range:

- Remains below the T_g
- Maintains chain stability
- Does not alter surface tension

CONSEQUENTLY:

- Micro-crack formation is prevented
- Optical degradation does not occur
- Internal lumen surfaces remain stable

POLYMER DAMAGE: CLINICAL SIGNIFICANCE

Micro-Damage Mechanism

High temperature and oxidative stress lead to:

- Capillary crack formation
- Surface roughening
- Protein adhesion

Clinical Effects

Such damage results in:

- Biofilm formation
- Sterilization resistance
- Increased infection risk

HDROZONE™ APPROACH



The HDROZONE™ approach operates on the basis of:

- Low temperature
- Controlled oxidation
- Stabilized reactive species

Clinical effects:

- Surface is preserved
- Device service life is extended
- Re-sterilization safety is enhanced

STERILIZATION WITHIN THE LUMEN



Critical Clinical Problem

The most vulnerable component of endoscopes is

- The internal lumen surface.

Limitations of Current Systems

- Efficacy diminishes as gas advances through the lumen
- Sterility decreases at the distal end


HDROZONE Approach

The HDROZONE™ approach achieves intraluminal performance through:

- Directed delivery of reactive species
- Central transport within the lumen
- Reduced surface loss

RESULT:

- Deep penetration
- Homogeneous sterilization



Residuals and Toxicity

Post-sterilization residuals can cause:

- Tissue irritation
- Toxic reactions
- Toxic Anterior Segment Syndrome (TASS)




The HDROZONE™ approach addresses this through:

- Decomposition of reactive species at the end of the process
- Reduction of residuals to a minimum level

RESULT:

- Patient safety is improved
- Risk of complications is reduced

COMPARATIVE CLINICAL ASSESSMENT



Key market participants in low-temperature sterilization include Advanced Sterilization Products (ASP) and STERIS.

| Criterion | Conventional H ₂ O ₂ Systems | HdrOzone |
|-------------------|--|------------|
| Temperature | 45–55 °C | 37–42 °C |
| Polymer impact | Potentially | Controlled |
| Lumen performance | Damaging | Enhanced |
| Residuals | Limited | Minimized |
| Humidity effect | May be present | Tolerable |

CLINICAL USE BENEFITS

Endoscope Use

- Extended device service life
- More stable image quality
- Safe reuse

Surgical Safety

- Reduced infection risk
- Enhanced sterilization reliability

Operational Continuity

- Fewer device failures
- Reduced need for re-sterilization

IMPACT ON CLINICAL OUTCOMES

The HDROZONE™ approach:

- Improves patient safety
- Reduces complications
- Enhances device performance

FUTURE PERSPECTIVE

Sterilization technologies are evolving away from high-temperature-based systems toward controlled-chemistry, low-temperature-based systems.

CONCLUSION

HDROZONE™:

- Operates at ultra-low temperature
- Protects polymers
- Enhances lumen performance
- Minimizes residuals

OFFICIAL CONCLUSION STATEMENT

This study was prepared with the objective of assessing whether the HDROZONE™ sterilization approach can achieve microbiological efficacy in flexible endoscopes while remaining within acceptable limits with respect to material integrity, optical quality, functional performance, and residual safety, and to establish the technical framework for defining the non-damaging sterilization process window at ultra-low temperature.

REFERENCES

Bu Daniela Trogolo, J. Samuel Arey, Peter R. Tentscher. Gas-Phase Ozone Reactions with a Structurally Diverse Set of Molecules: Barrier Heights and Reaction Energies Evaluated by Coupled Cluster and Density Functional Theory Calculations. *The Journal of Physical Chemistry A* 2019, 123 (2) , 517-536. DOI: 10.1021/acs.jpca.8b10323.

1. Alexandra Fischbacher, Katja Löppenber, Clemens von Sonntag, and Torsten C. Schmidt . A New Reaction Pathway for Bromite to Bromate in the Ozonation of Bromide. *Environmental Science & Technology* 2015, 49 (19) , 11714-11720. DOI: 10.1021/acs.est.5b02634.

2. Yi Yang, Jin Jiang, Jinglin Lu, Jun Ma, and Yongze Liu . Production of Sulfate Radical and Hydroxyl Radical by Reaction of Ozone with Peroxymonosulfate: A Novel Advanced Oxidation Process. *Environmental Science & Technology* 2015, 49 (12) , 7330-7339. DOI: 10.1021/es506362e.

3. Janez Cerkovnik and Božo Plesničar . Recent Advances in the Chemistry of Hydrogen Trioxide (HOOH). *Chemical Reviews* 2013, 113 (10) , 7930-7951. DOI: 10.1021/cr300512s.

4. Alexandra Fischbacher, Justus von Sonntag, Clemens von Sonntag, and Torsten C. Schmidt . The •OH Radical Yield in the H₂O₂ + O₃ (Peroxone) Reaction. *Environmental Science & Technology* 2013, 47 (17) , 9959-9964. DOI: 10.1021/es402305r.

5. Carlos Barrera-Díaz, Lina. A. Bernal-Martínez, Reyna Natividad, and Juan M. Peralta-Hernández . Synergy of Electrochemical/O₃ Process with Aluminum Electrodes in Industrial Wastewater Treatment. *Industrial & Engineering Chemistry Research* 2012, 51 (27) , 9335-9342. DOI: 10.1021/ie3004144.

6. Sergej Naumov and Clemens von Sonntag . Standard Gibbs Free Energies of Reactions of Ozone with Free Radicals in Aqueous Solution: Quantum-Chemical Calculations. *Environmental Science & Technology* 2011, 45 (21) , 9195-9204. DOI: 10.1021/es2018658.

7. J. Pablo Pocostales, Myint M. Sein, Wolfgang Knolle, Clemens von Sonntag, and Torsten C. Schmidt. Degradation of Ozone-Refractory Organic Phosphates in Wastewater by Ozone and Ozone/Hydrogen Peroxide (Peroxone): The Role of Ozone Consumption by Dissolved Organic Matter. *Environmental Science & Technology* 2010, 44 (21) , 8248-8253. DOI: 10.1021/es1018288.

8. Helena Jablonowski, Joao Santos Sousa, Klaus-Dieter Weltmann, Kristian Wende, Stephan Reuter. Quantification of the ozone and singlet delta oxygen produced in gas and liquid phases by a non-thermal atmospheric plasma with relevance for medical treatment. *Scientific Reports* 2018, 8 (1) DOI: 10.1038/s41598-018-30483-w.

9. Yalei Ding, Jiejie Wang, Shanshan. u, Kun-Yi Andrew Lin, Shaoping Tong. Oxygen vacancy of CeO₂ improved efficiency of H₂O₂/O₃ for the degradation of acetic acid in acidic solutions. *Separation and Purification Technology* 2018, 207, 92-98. DOI: 10.1016/j.seppur.2018.06.027.

10. Donggwan Lee, Jae-Cheol Lee, Joo-Youn Nam, Hyun-Woo Kim. Degradation of sulfonamide antibiotics and their intermediates toxicity in an aeration-assisted non-thermal plasma while treating strong wastewater. *Chemosphere* 2018, 209, 901-907. DOI: 10.1016/j.chemosphere.2018.06.125.
11. David B. Miklos, Christian Remy, Martin Jekel, Karl G. Linden, Jörg E. Drewes, Uwe Hübner. Evaluation of advanced oxidation processes for water and wastewater treatment – A critical review. *Water Research* 2018, 139, 118-131. DOI: 10.1016/j.watres.2018.03.042.
12. Pierre-François Biard, Thom Thi Dang, Jenny Bocanegra, Annabelle Couvert. Intensification of the O₃/H₂O₂ advanced oxidation process using a continuous tubular reactor filled with static mixers: Proof of concept. *Chemical Engineering Journal* 2018, 344, 574-582. DOI: 10.1016/j.cej.2018.03.112.
13. Emmanuel Mousset, Nihal Oturan, Mehmet A. Oturan. An unprecedented route of OH radical reactivity evidenced by an electrocatalytical process: Ipso-substitution with perhalogenocarbon compounds. *Applied Catalysis B: Environmental* 2018, 226, 135-146. DOI: 10.1016/j.apcatb.2017.12.028.
14. J.F. Pérez, S. Sabatino, A. Galia, M.A. Rodrigo, J. Llanos, C. Sáez, O. Scialdone. Effect of air pressure on the electro-Fenton process at carbon felt electrodes. *Electrochimica Acta* 2018, 273, 447-453. DOI: 10.1016/j.electacta.2018.04.031.
15. Guilherme Garcia Bessegato, João Carlos de Souza, Juliano Carvalho Cardoso, Maria Valnice Boldrin Zanoni. Assessment of several advanced oxidation processes applied in the treatment of environmental concern constituents from a real hair dye wastewater. *Journal of Environmental Chemical Engineering* 2018, 6 (2) , 2794-2802. DOI: 10.1016/j.jece.2018.04.041.
16. A. Privat-Maldonado, Y. Gorbanev, D. O'Connell, R. Vann, V. Chechik, M. W. van der Woude. Nontarget Biomolecules Alter Macromolecular Changes Induced by Bactericidal Low-Temperature Plasma. *IEEE Transactions on Radiation and Plasma Medical Sciences* 2018, 2 (2) , 121-128. DOI: 10.1109/TRPMS.2017.2761405.
17. Mika Sillanpää, Mohamed Chaker Ncibi, Anu Matilainen. Advanced oxidation processes for the removal of natural organic matter from drinking water sources: A comprehensive review. *Journal of Environmental Management* 2018, 208, 56-76. DOI: 10.1016/j.jenvman.2017.12.009.
18. B.W. Darvell. *More Chemistry*. 2018,,, 771-789. DOI: 10.1016/B978-0-08-101035-8.50030-4.
19. Erika Reisz, Clemens von Sonntag, Agnes Tekle-Röttering, Sergej Naumov, Winfried Schmidt, Torsten C. Schmidt. Reaction of 2-propanol with ozone in aqueous media. *Water Research* 2018, 128, 171-182. DOI: 10.1016/j.watres.2017.10.035.
20. João F. Gomes, Inês Leal, Katarzyna Bednarczyk, Marta Gmurek, Marek Stelmachowski, Magdalena Diak, M. Emília Quinta-Ferreira, Raquel Costa, Rosa M. Quinta-Ferreira, Rui C. Martins. Photocatalytic ozonation using doped TiO₂ catalysts for the removal of parabens in water. *Science of The Total Environment* 2017, 609, 329-340. DOI: 10.1016/j.scitotenv.2017.07.180.

21. Pierre-François Biard, Thom Thi Dang, Annabelle Couvert. Determination by reactive absorption of the rate constant of the ozone reaction with the hydroperoxide anion. *Chemical Engineering Research and Design* 2017, 127, 62-71. DOI: 10.1016/j.cherd.2017.09.004.
22. Ozge Turkay, Sibel Barışçı, Mika Sillanpää. E-peroxone process for the treatment of laundry wastewater: A case study. *Journal of Environmental Chemical Engineering* 2017, 5 (5) , 4282-4290. DOI: 10.1016/j.jece.2017.08.012.
23. Minhwan Kwon, Homin Kye, Youmi Jung, Yejoon Yoon, Joon-Wun Kang. Performance characterization and kinetic modeling of ozonation using a new method: R OH,O₃ concept. *Water Research* 2017, 122, 172-182. DOI: 10.1016/j.watres.2017.05.062.
24. Ligy Philip, Bhallamudi Murty, Channarong Puchongkawarin, Miao Guo, Nilay Shah, David Stuckey, Benoit Chachuat, Yannic Vaupel, Sarojini Tiwari, Chitta Behera, Babji Srinivasan, Chedly Tizaoui, Olajumoke Odejimi, Ayman Abdelaziz. *3 Wastewater Treatment*. 2017,, 115-180. DOI: 10.1201/9781315153209-4.
25. Songjie Wu, Qian Zhang, Ruonan Ma, Shuang Yu, Kaile Wang, Jue Zhang, Jing Fang. Reactive radical-driven bacterial inactivation by hydrogen-peroxide-enhanced plasma-activated-water. *The European Physical Journal Special Topics* 2017, 226 (13) , 2887-2899. DOI: 10.1140/epjst/e2016-60330-y.
26. Julia Patzsch, Jonathan Z. Bloh. Improved photocatalytic ozone abatement over transition metal-grafted titanium dioxide. *Catalysis Today* 2017, DOI: 10.1016/j.cattod.2017.07.010.
27. Xinyang Li, Shaobin Sun., u Zhang, Guicheng Liu, Clark Renjun Zheng, Jianzhong Zheng, Dayi Zhang, Hong Yao. Combined electro-catazone/electro-peroxone process for rapid and effective Rhodamine B degradation. *Separation and Purification Technology* 2017, 178, 189-192. DOI: 10.1016/j.seppur.2016.12.052.
28. Jingxin Yang, Ji Li, Wenyi Dong, Jun Ma, Jiayin Li. Influence of nitrite on the degradation of atrazine by ozonation. *Journal of Chemical Technology & Biotechnology* 2017, 92 (2) , 442-450. DOI: 10.1002/jctb.5031.
29. George Wafula Wanjala, Arnold Onyango, Calvin Onyango, Moses Makayoto. . *African Journal of Biochemistry Research* 2017,, 79. DOI: 10.5897/AJBR2017.0967.
30. Glen Andrew de Vera, Wolfgang Gernjak, Howard Weinberg, Maria José Farré, Jurg Keller, Urs von Gunten. Kinetics and mechanisms of nitrate and ammonium formation during ozonation of dissolved organic nitrogen. *Water Research* 2017, 108, 451-461. DOI: 10.1016/j.watres.2016.10.021.
31. M. E. Zappi, R. Hernandez, D. Gang, R. Bajpai, C. H. Kuo, D. O. Hill. Treatment of groundwater contaminated with high levels of explosives using advanced oxidation processes. *International Journal of Environmental Science and Technology* 2016, 13 (12) , 2767-2778. DOI: 10.1007/s13762-016-1109-x.
32. Juhong Zhan, Yujue Wang, Huijiao Wang, Wenhua Shen., uejun Pan, Jinlin Wang, Gang Yu. Electro-peroxone regeneration of phenol-saturated activated carbon fiber: The effects of irreversible adsorption and operational parameters. *Carbon* 2016, 109, 321-330. DOI: 10.1016/j.carbon.2016.08.034.

33. Jingxin Yang, Ji Li, Wenyi Dong, Jun Ma, Jie Cao, Tingting Li, Jiayin Li, Jia Gu, Pingxin Liu. Study on enhanced degradation of atrazine by ozonation in the presence of hydroxylamine. *Journal of Hazardous Materials* 2016, 316, 110-121. DOI: 10.1016/j.jhazmat.2016.04.078.
34. Xin Cheng, Hongguang Guo, Hongwei Liu, Yang Liu, Ying Yang, Yongli Zhang. Performance and Mechanism on Degradation of Estriol Using O₃/PS Process. *Ozone: Science & Engineering* 2016, 38 (5) , 358-366. DOI: 10.1080/01919512.2016.1170589.
35. László Wojnárovits, Erzsébet Takács. Radiation Induced Degradation of Organic Pollutants in Waters and Wastewaters. *Topics in Current Chemistry* 2016, 374 (4) DOI: 10.1007/s41061-016-0050-2.
36. Jaedon Shin, Zahra Ramadhany Hidayat, Yunho Lee. Influence of Seasonal Variation of Water Temperature and Dissolved Organic Matter on Ozone and OH Radical Reaction Kinetics During Ozonation of a Lake Water. *Ozone: Science & Engineering* 2016, 38, 100-114. DOI: 10.1080/01919512.2015.1079120.
37. Yury Gorbanev, Deborah O'Connell, Victor Chechik. Non-Thermal Plasma in Contact with Water: The Origin of Species. *Chemistry - A European Journal* 2016, 22 (10) , 3496-3505. DOI: 10.1002/chem.201503771.
38. Arnold N. Onyango. Alternatives to the 'water oxidation pathway' of biological ozone formation. *Journal of Chemical Biology* 2016, 9, 1-8. DOI: 10.1007/s12154-015-0140-6.
39. Fei Qi, Wei Chu, Bingbing. u. Comparison of phenacetin degradation in aqueous solutions by catalytic ozonation with CuFe₂O₄ and its precursor: Surface properties, intermediates and reaction mechanisms. *Chemical Engineering Journal* 2016, 284, 28-36. DOI: 10.1016/j.cej.2015.07.095.
40. Erika Reisz, Sergej Naumov, Winfried Schmidt, Clemens von Sonntag. Reaction of Ozone with Ag(I)—Mechanistic Considerations. *Ozone: Science & Engineering* 2015, 37, 393-404. DOI: 10.1080/01919512.2015.1041583.
41. Huijiao Wang, Shi Yuan, Juhong Zhan, Yujue Wang, Gang Yu, Shubo Deng, Jun Huang, Bin Wang. Mechanisms of enhanced total organic carbon elimination from oxalic acid solutions by electro-peroxone process. *Water Research* 2015, 80, 20-29. DOI: 10.1016/j.watres.2015.05.024.
42. Jesús Ferre-Aracil, Salvador C. Cardona, Javier Navarro-Laboulais. Kinetic Study of Ozone Decay in Homogeneous Phosphate-Buffered Medium. *Ozone: Science & Engineering* 2015, 37, 330-342. DOI: 10.1080/01919512.2014.998756.
43. Fei Qi, Bingbing. u, Wei Chu. Heterogeneous catalytic ozonation of phenacetin in water using magnetic spinel ferrite as catalyst: Comparison of surface property and efficiency. *Journal of Molecular Catalysis A: Chemical* 2015, 396, 164-173. DOI: 10.1016/j.molcata.2014.10.001.
44. Erika Reisz, Alexandra Fischbacher, Sergej Naumov, Clemens von Sonntag, Torsten C. Schmidt. Hydride Transfer: A Dominating Reaction of Ozone with Tertiary Butanol and Formate Ion in Aqueous Solution. *Ozone: Science & Engineering* 2014, 36, 532-539. DOI: 10.1080/01919512.2014.891436.
45. Sergey L. Khursan. Peroxide intermediates of oxidation processes: Organic trioxides. 2014,, 1-72. DOI: 10.1002/9780470682531.pat0874.

46. S. Reed Plimpton, Mark Gołkowski, Deborah G. Mitchell, Chad Austin, Sandra S. Eaton, Gareth R. Eaton, Czesław Gołkowski, Martin Voskuil. Remote delivery of hydroxyl radicals via secondary chemistry of a nonthermal plasma effluent. *Biotechnology and Bioengineering* 2013, 110 (10.1002/bit.v110.7) , 1936-1944. DOI: 10.1002/bit.24853.
47. Zdenko Machala, Barbora Tarabova, Karol Hensel, Eva Spetlikova, Libusa Sikurova, Petr Lukes. Formation of ROS and RNS in Water Electro-Sprayed through Transient Spark Discharge in Air and their Bactericidal Effects. *Plasma Processes and Polymers* 2013, 10 (10.1002/ppap.v10.7) , 649-659. DOI: 10.1002/ppap.201200113.
48. Josep M. Anglada, Miquel Torrent-Sucarrat, Manuel F. Ruiz-Lopez, Marilia Martins-Costa. Is the HO₄⁻ Anion a Key Species in the Aqueous-Phase Decomposition of Ozone?. *Chemistry - A European Journal* 2012, 18 (10.1002/chem.v18.42) , 13435-13445. DOI: 10.1002/chem.201201991.
49. Marek Golkowski, Czesław Golkowski, Jori Leszczynski, S. Reed Plimpton, Piotr Maslowski, Aleksandra Foltynowicz, Jun Ye, Bruce McCollister. Hydrogen-Peroxide-Enhanced Nonthermal Plasma Effluent for Biomedical Applications. *IEEE Transactions on Plasma Science* 2012, 40, 1984-1991. DOI: 10.1109/TPS.2012.2200910.
50. D. Gardoni, A. Vailati, R. Canziani. Decay of Ozone in Water: A Review. *Ozone: Science & Engineering* 2012, 34, 233-242. DOI: 10.1080/01919512.2012.686354.
51. Sergej Naumov, Clemens von Sonntag. The reaction of •OH with O₂, the decay of O₃^{•-} and the pK_a of HO₃^{•-} - interrelated questions in aqueous free-radical chemistry. *Journal of Physical Organic Chemistry* 2011, 24 (10.1002/poc.v24.7) , 600-602. DOI: 10.1002/poc.1812.
52. Rossano Amadelli, Luca Samiolo, Achille De Battisti, Alexander B. Velichenko. Electro-oxidation of Some Phenolic Compounds by Electrogenenerated O₃ and by Direct Electrolysis at PbO₂ Anodes. *Journal of The Electrochemical Society* 2011, 158, P87. DOI: 10.1149/1.3589913.
53. Pentax, White Papers
54. ASP, White Papers
55. Olympus, White Papers