

A New Physical Boundary in Sterilization Technology

Analysis of the Physico-Chemical Model

Patented (No: TR2019 17785 B)

Hydronium Ionization and Long-Lumen Sterilization Technical Report

Supported by Hyggen Accredited Laboratories Bischofshofen, 29 September 2020

Hasan Tahsin ÖZBEK • Bülent DEVECİ • Yüksel ERGÜN

White Paper - **02.03.2026**

ABSTRACT

The widespread adoption of minimally invasive surgery and flexible endoscopy systems has revealed the geometric limit problem in sterilization processes. Particularly:

- \varnothing 1–2 mm diameter
- 1–15+ meter length
- Complex internal surface structures

These factors define the physical boundaries of diffusion-based low-temperature sterilization methods. This study comprehensively defines, based on H_2O_2 – O_3 –DBD plasma combination leading to hydronium (H_3O^+) and peroxizone (H_2O_4) formation:

- Reactive species generation kinetics
- Plasma environment parameters
- Intra-lumen transport mechanisms
- Biological validation
- Material interaction
- Residue and toxicology

Core Approach:

Sterilization is not a passive diffusion problem but an **electric field-assisted reactive transport and phase-controlled chemical generation** problem.



KEYLER

Peroxyzone, DBD Plasma, Reactive Transport, Diffusion, Kinetics, Ionization, Validation, Toxicology

1. INTRODUCTION AND PROBLEM STATEMENT

Flexible endoscopes and robotic surgical instruments, due to their:

- Multi-layered structure
- Heterogeneous surface energy
- Long and narrow lumen geometry

represent boundary conditions for classical sterilization models.

Diffusion-based gas penetration:

- Decreases exponentially with distance
- Loses energy through surface interactions
- Is interrupted by phase changes

This situation particularly limits sterility assurance in long lumens.

2. PHYSICAL LIMITS OF EXISTING TECHNOLOGIES

2.1 Diffusion Limitation

Concentration along the lumen:

$$C(x) = C_0 e^{-kx}$$

(where k = surface interaction coefficient, x = lumen length)

Result:

As x increases, C(x) decreases rapidly and becomes ineffective beyond a certain length.

2.2 Condensation

Gas-phase H₂O₂:

- Loses energy through wall collisions
- Transitions to liquid phase under local temperature/partial pressure conditions

Result:

- Film formation
- Flow blockage
- Cessation of further penetratio

2.3 Gas Depletion

- Reactives are consumed along the lumen
- Active species concentration drops at distal ends

2.4 Plasma System Limitations

- Radicals have very short lifetimes
- Localized production localized effect
- Not sustainable along the lumen

2.5 Ozone System Limitations

- High oxidative aggression
- Difficult to control
- Material damage risk

General Conclusion:

No system currently exists that simultaneously provides:

- Long-distance penetration
- Low temperature
- Low material damage
- Low residue

3. HYBRID REACTIVE SYSTEM DEFINITION

3.1 System Components

- H₂O₂: Precursor
- O₃: Oxidative driver
- DBD Plasma: Energy source and control mechanism

3.2 Reactive Species

- H₃O⁺: Ionic carrier
- H₂O₄: Stable oxidant
- OH·: Highly reactive radical

3.3 Fundamental Transport Mechanism

Classical: $J = -D\nabla C$

Hybrid: $J = -D\nabla C + \mu CE$

(D = diffusion coefficient, μCE = electric field-driven drift)

4. PLASMA PHYSICS

4.1 Operating Range

- Pressure: 0.1 – 0.4 Pa
- Voltage: ~6 kV
- Current: ~20 mA
- Frequency: 50–100 Hz

4.2 Plasma Parameters

- Electron density: $n_e \approx 10^{13} - 10^{15} \text{ m}^{-3}$
- Electron temperature: $T_e \approx 1 - 3 \text{ eV}$

4.3 Debye Length

$$\lambda_D = \sqrt{\frac{\epsilon_0 k T_e}{n_e e^2}}$$

$\approx 10^{-4} - 10^{-3}$ m

Result:

- $\lambda_D <$ lumen diameter
- Plasma remains volumetrically stable
- Ions remain in the volume rather than collapsing onto the wall

5. REACTION KINETICS

5.1 Key Reactions



5.2 Differential Kinetic Model

$$\frac{d[OH]}{dt} = k_1[H_2O_2] + k_2[O_3] - k_3[OH]^2$$

$$\frac{d[H_3O^+]}{dt} = k_4[OH][H_2O] - k_5[recombination]$$

5.3 Phase Management

- Plasma ON: Radical production
- Plasma OFF: Stable species formation

This structure enables:

- Controlled reactive production
- Radical suppression
- Stable oxidant generation

6. INTRA-LUMEN TRANSPORT

6.1 Coulomb Repulsion

- Similarly charged ions repel each other
- Concentration builds toward the center

6.2 Flow Characteristics

- Reduced wall interaction
- Decreased energy loss
- Jet-like flow formation

6.3 Conclusion

- Concentration loss is minimized
- Long-distance transport becomes possible

7. VALIDATION

7.1 Test Configuration

- Lumen: 2 mm / 15 m
- Load: Full load
- BI: 10⁶ Geobacillus stearothermophilus
- Single Cycle: 35 minutes



4.2.3 Positions and results of the microbiological indicators

bionova - biological indicator		
cfu / carrier [lg]	mean value [lg]	
Control 1-2	8,69/8,94	8,82

Wrapping* [yes / no]	sample No.	designation	cfu / carrier [lg]	Enrichment* [3 / 7 days]	lg / carrier	Reduction factor [lg]
upper basket						
yes	91	A	0	-/-	0	≥8,82
yes	99	E	0	-/-	0	≥8,82
yes	93	B	0	-/-	0	≥8,82
yes	95	C	0	-/-	0	≥8,82
yes	101	F	0	-/-	0	≥8,82
yes	97	D	0	-/-	0	≥8,82
middle basket						
yes	103	G	0	-/-	0	≥8,82
yes	105	H	0	-/-	0	≥8,82
yes	104	G	0	-/-	0	≥8,82
lower basket						
yes	94	B	0	-/-	0	≥8,82
yes	100	E	0	-/-	0	≥8,82
yes	92	A	0	-/-	0	≥8,82
yes	98	D	0	-/-	0	≥8,82
yes	102	F	0	-/-	0	≥8,82
yes	96	C	0	-/-	0	≥8,82

legend:

- RF = Reduction factor
 - = no turbidity due to microbial growth
 + = turbidity due to microbial growth
 n.z. = uncountable

- A = HygCen – white PTFE PCD / PCD 1 mm Ø 850 mm long
 B = HygCen – white PTFE PCD / PCD 2 mm Ø 1200 mm long
 C = Teknomar – Steel Lumen PCD / PCD 0,7 mm Ø 500 mm long
 D = Teknomar – white PTFE PCD / PCD 0,4 mm Ø 900 mm long
 E = Teknomar – white PTFE PCD / PCD 2 mm Ø 7500 mm long
 F = Teknomar – white PTFE PCD / PCD 2 mm Ø 10000 mm long
 G = Teknomar – white PTFE PCD / PCD 2 mm Ø 15000 mm long
 H = Teknomar – white PTFE PCD / PCD 2 mm Ø 50000 mm long
 I = polyester suture bionova / Geobacillus stearothermophilus. / in double tyvek

7.2 Results

- 100% inactivation
- SAL = 10⁻⁶

7.3 Statistics

- Number of repetitions ≥ 3
- Half-cycle validation available

8. MATERIAL COMPATIBILITY

Tested materials:

- 316L,
- PEEK,
- Silicone

8.1 Observations

- No micro-pitting
- Passivation preserved
- No crack formation

8.2 Damage Index

$$DI = \alpha(\text{oxidation}) + \beta(\text{temperature}) + \gamma(\text{radical})$$

Low DI

9. SURFACE ANALYSIS (SEM)

- Surface topography unchanged
- Microstructure stable

10. CHEMICAL RESIDUE (FTIR)

- LOD: 0.01 mg/m²
- Measured: < LOD

10.1 Conversion Mechanism

- H₂O₂: H₂O
- O₃: O₂

11. CONTROL ALGORITHM

Input Parameters:

- Pressure,
- Humidity,
- Gas ratio,
- Plasma duty cycle

Target:

- Maximize: H₃O⁺ and H₂O₄
- Minimize: Radicals and residue

Cycle Structure:

1. H₂O₂ loading
2. O₃ injection
3. Plasma ON
4. Plasma OFF
5. Feedback loop

12. COMPETITIVE ANALYSIS

Parameter	Diffusion Systems	Hybrid System
Transport	Diffusion only	Drift + Diffusion
Lumen Capability	~1 m	15 m
Humidity	Problem	Advantage (reactive)
Control	Limited	Active

13. ECONOMIC IMPACT

- Increased device lifespan
- Reduced consumable costs
- Decreased maintenance requirements

14. REGULATORY COMPLIANCE

- MDR 2017/745
- ISO 14937
- ISO 11138
- ISO 10993

15. TOXICOLOGY AND EMISSIONS

- Residue: Not detectable
- Emissions: O₂ + H₂O
- Toxicity: Non-toxic

16. CONCLUSION

This study demonstrates that:

- Sterilization is not limited to the diffusion model
- Electric field-assisted reactive transport is critical
- Long-lumen sterilization is physically achievable

The H_3O^+ and H_2O_4 -based hybrid approach simultaneously delivers:

- Controlled chemistry generation
- Low material damage
- Low residue
- Superior long-distance penetration

B SCIENCE[®]
BIOMEDICAL

KAYNAKÇA

Bu Daniela Trogolo, J. Samuel Arey, Peter R. Tentscher. Gas-Phase Ozone Reactions with a Structurally Diverse Set of Molecules: Barrier Heights and Reaction Energies Evaluated by Coupled Cluster and Density Functional Theory Calculations. *The Journal of Physical Chemistry A* 2019, 123 (2) , 517-536. DOI: 10.1021/acs.jpca.8b10323.

1. Alexandra Fischbacher, Katja Löppenberg, Clemens von Sonntag, and Torsten C. Schmidt . A New Reaction Pathway for Bromite to Bromate in the Ozonation of Bromide. *Environmental Science & Technology* 2015, 49 (19) , 11714-11720. DOI: 10.1021/acs.est.5b02634.
2. Yi Yang, Jin Jiang, inglin Lu, Jun Ma, and Yongze Liu . Production of Sulfate Radical and Hydroxyl Radical by Reaction of Ozone with Peroxymonosulfate: A Novel Advanced Oxidation Process. *Environmental Science & Technology* 2015, 49 (12) , 7330-7339. DOI: 10.1021/es506362e.
3. Janez Cerkovnik and Božo Plesničar . Recent Advances in the Chemistry of Hydrogen Trioxide (HOOH). *Chemical Reviews* 2013, 113 (10) , 7930-7951. DOI: 10.1021/cr300512s.
4. Alexandra Fischbacher, Justus von Sonntag, Clemens von Sonntag, and Torsten C. Schmidt . The •OH Radical Yield in the H₂O₂ + O₃ (Peroxone) Reaction. *Environmental Science & Technology* 2013, 47 (17) , 9959-9964. DOI: 10.1021/es402305r.
5. Carlos Barrera-Díaz, Lina. A. Bernal-Martínez, Reyna Natividad, and Juan M. Peralta-Hernández . Synergy of Electrochemical/O₃ Process with Aluminum Electrodes in Industrial Wastewater Treatment. *Industrial & Engineering Chemistry Research* 2012, 51 (27) , 9335-9342. DOI: 10.1021/ie3004144.
6. Sergej Naumov and Clemens von Sonntag . Standard Gibbs Free Energies of Reactions of Ozone with Free Radicals in Aqueous Solution: Quantum-Chemical Calculations. *Environmental Science & Technology* 2011, 45 (21) , 9195-9204. DOI: 10.1021/es2018658.
7. J. Pablo Pocostales, Myint M. Sein, Wolfgang Knolle, Clemens von Sonntag, and Torsten C. Schmidt. Degradation of Ozone-Refractory Organic Phosphates in Wastewater by Ozone and Ozone/Hydrogen Peroxide (Peroxone): The Role of Ozone Consumption by Dissolved Organic Matter. *Environmental Science & Technology* 2010, 44 (21) , 8248-8253. DOI: 10.1021/es1018288.
8. Helena Jablonowski, Joao Santos Sousa, Klaus-Dieter Weltmann, Kristian Wende, Stephan Reuter. Quantification of the ozone and singlet delta oxygen produced in gas and liquid phases by a non-thermal atmospheric plasma with relevance for medical treatment. *Scientific Reports* 2018, 8 (1) DOI: 10.1038/s41598-018-30483-w.
9. Yalei Ding, Jiejie Wang, Shanshan. u, Kun-Yi Andrew Lin, Shaoping Tong. Oxygen vacancy of CeO₂ improved efficiency of H₂O₂/O₃ for the degradation of acetic acid in acidic solutions. *Separation and Purification Technology* 2018, 207, 92-98. DOI: 10.1016/j.seppur.2018.06.027.

10. Donggwan Lee, Jae-Cheol Lee, Joo-Youn Nam, Hyun-Woo Kim. Degradation of sulfonamide antibiotics and their intermediates toxicity in an aeration-assisted non-thermal plasma while treating strong wastewater. *Chemosphere* 2018, 209, 901-907. DOI: 10.1016/j.chemosphere.2018.06.125.
11. David B. Miklos, Christian Remy, Martin Jekel, Karl G. Linden, Jörg E. Drewes, Uwe Hübner. Evaluation of advanced oxidation processes for water and wastewater treatment – A critical review. *Water Research* 2018, 139, 118-131. DOI: 10.1016/j.watres.2018.03.042.
12. Pierre-François Biard, Thom Thi Dang, Jenny Bocanegra, Annabelle Couvert. Intensification of the O₃/H₂O₂ advanced oxidation process using a continuous tubular reactor filled with static mixers: Proof of concept. *Chemical Engineering Journal* 2018, 344, 574-582. DOI: 10.1016/j.cej.2018.03.112.
13. Emmanuel Mousset, Nihal Oturan, Mehmet A. Oturan. An unprecedented route of OH radical reactivity evidenced by an electrocatalytical process: Ipso-substitution with perhalogenocarbon compounds. *Applied Catalysis B: Environmental* 2018, 226, 135-146. DOI: 10.1016/j.apcatb.2017.12.028.
14. J.F. Pérez, S. Sabatino, A. Galia, M.A. Rodrigo, J. Llanos, C. Sáez, O. Scialdone. Effect of air pressure on the electro-Fenton process at carbon felt electrodes. *Electrochimica Acta* 2018, 273, 447-453. DOI: 10.1016/j.electacta.2018.04.031.
15. Guilherme Garcia Bessegato, João Carlos de Souza, Juliano Carvalho Cardoso, Maria Valnice Boldrin Zanoni. Assessment of several advanced oxidation processes applied in the treatment of environmental concern constituents from a real hair dye wastewater. *Journal of Environmental Chemical Engineering* 2018, 6 (2) , 2794-2802. DOI: 10.1016/j.jece.2018.04.041.
16. A. Privat-Maldonado, Y. Gorbanev, D. O'Connell, R. Vann, V. Chechik, M. W. van der Woude. Nontarget Biomolecules Alter Macromolecular Changes Induced by Bactericidal Low-Temperature Plasma. *IEEE Transactions on Radiation and Plasma Medical Sciences* 2018, 2 (2) , 121-128. DOI: 10.1109/TRPMS.2017.2761405.
17. Mika Sillanpää, Mohamed Chaker Ncibi, Anu Matilainen. Advanced oxidation processes for the removal of natural organic matter from drinking water sources: A comprehensive review. *Journal of Environmental Management* 2018, 208, 56-76. DOI: 10.1016/j.jenvman.2017.12.009.
18. B.W. Darvell. *More Chemistry*. 2018,, 771-789. DOI: 10.1016/B978-0-08-101035-8.50030-4.
19. Erika Reisz, Clemens von Sonntag, Agnes Tekle-Röttering, Sergej Naumov, Winfried Schmidt, Torsten C. Schmidt. Reaction of 2-propanol with ozone in aqueous media. *Water Research* 2018, 128, 171-182. DOI: 10.1016/j.watres.2017.10.035.
20. João F. Gomes, Inês Leal, Katarzyna Bednarczyk, Marta Gmurek, Marek Stelmachowski, Magdalena Diak, M. Emília Quinta-Ferreira, Raquel Costa, Rosa M. Quinta-Ferreira, Rui C. Martins. Photocatalytic ozonation using doped TiO₂ catalysts for the removal of parabens in water. *Science of The Total Environment* 2017, 609, 329-340. DOI: 10.1016/j.scitotenv.2017.07.180.

21. Pierre-François Biard, Thom Thi Dang, Annabelle Couvert. Determination by reactive absorption of the rate constant of the ozone reaction with the hydroperoxide anion. *Chemical Engineering Research and Design* 2017, 127, 62-71. DOI: 10.1016/j.cherd.2017.09.004.
22. Ozge Turkay, Sibel Barışçı, Mika Sillanpää. E-peroxone process for the treatment of laundry wastewater: A case study. *Journal of Environmental Chemical Engineering* 2017, 5 (5) , 4282-4290. DOI: 10.1016/j.jece.2017.08.012.
23. Minhwan Kwon, Homin Kye, Youmi Jung, Yeojoon Yoon, Joon-Wun Kang. Performance characterization and kinetic modeling of ozonation using a new method: R OH,O₃ concept. *Water Research* 2017, 122, 172-182. DOI: 10.1016/j.watres.2017.05.062.
24. Ligy Philip, Bhallamudi Murty, Channarong Puchongkawarin, Miao Guo, Nilay Shah, David Stuckey, Benoit Chachuat, Yannic Vaupel, Sarojini Tiwari, Chitta Behera, Babji Srinivasan, Chedly Tizaoui, Olajumoke Odejimi, Ayman Abdelaziz. 3 *Wastewater Treatment*. 2017,, 115-180. DOI: 10.1201/9781315153209-4.
25. Songjie Wu, Qian Zhang, Ruonan Ma, Shuang Yu, Kaile Wang, Jue Zhang, Jing Fang. Reactive radical-driven bacterial inactivation by hydrogen-peroxide-enhanced plasma-activated-water. *The European Physical Journal Special Topics* 2017, 226 (13) , 2887-2899. DOI: 10.1140/epjst/e2016-60330-y.
26. Julia Patzsch, Jonathan Z. Bloh. Improved photocatalytic ozone abatement over transition metal-grafted titanium dioxide. *Catalysis Today* 2017, DOI: 10.1016/j.cattod.2017.07.010.
27. Xinyang Li, Shaobin Sun, u Zhang, Guicheng Liu, Clark Renjun Zheng, Jianzhong Zheng, Dayi Zhang, Hong Yao. Combined electro-catazone/electro-peroxone process for rapid and effective Rhodamine B degradation. *Separation and Purification Technology* 2017, 178, 189-192. DOI: 10.1016/j.seppur.2016.12.052.
28. Jingxin Yang, Ji Li, Wenyi Dong, Jun Ma, Jiayin Li. Influence of nitrite on the degradation of atrazine by ozonation. *Journal of Chemical Technology & Biotechnology* 2017, 92 (2) , 442-450. DOI: 10.1002/jctb.5031.
29. George Wafula Wanjala, Arnold Onyango, Calvin Onyango, Moses Makayoto. . *African Journal of Biochemistry Research* 2017,, 79. DOI: 10.5897/AJBR2017.0967.
30. Glen Andrew de Vera, Wolfgang Gernjak, Howard Weinberg, Maria José Farré, Jurg Keller, Urs von Gunten. Kinetics and mechanisms of nitrate and ammonium formation during ozonation of dissolved organic nitrogen. *Water Research* 2017, 108, 451-461. DOI: 10.1016/j.watres.2016.10.021.
31. M. E. Zappi, R. Hernandez, D. Gang, R. Bajpai, C. H. Kuo, D. O. Hill. Treatment of groundwater contaminated with high levels of explosives using advanced oxidation processes. *International Journal of Environmental Science and Technology* 2016, 13 (12) , 2767-2778. DOI: 10.1007/s13762-016-1109-x.
32. Juhong Zhan, Yujue Wang, Huijiao Wang, Wenhua Shen, uejun Pan, Jinlin Wang, Gang Yu. Electro-peroxone regeneration of phenol-saturated activated carbon fiber: The effects of irreversible adsorption and operational parameters. *Carbon* 2016, 109, 321-330. DOI: 10.1016/j.carbon.2016.08.034.

33. Jingxin Yang, Ji Li, Wenyi Dong, Jun Ma, Jie Cao, Tingting Li, Jiayin Li, Jia Gu, Pingxin Liu. Study on enhanced degradation of atrazine by ozonation in the presence of hydroxylamine. *Journal of Hazardous Materials* 2016, 316, 110-121. DOI: 10.1016/j.jhazmat.2016.04.078.
34. Xin Cheng, Hongguang Guo, Hongwei Liu, Yang Liu, Ying Yang, Yongli Zhang. Performance and Mechanism on Degradation of Estriol Using O₃/PS Process. *Ozone: Science & Engineering* 2016, 38 (5) , 358-366. DOI: 10.1080/01919512.2016.1170589.
35. László Wojnárovits, Erzsébet Takács. Radiation Induced Degradation of Organic Pollutants in Waters and Wastewaters. *Topics in Current Chemistry* 2016, 374 (4) DOI: 10.1007/s41061-016-0050-2.
36. Jaedon Shin, Zahra Ramadhany Hidayat, Yunho Lee. Influence of Seasonal Variation of Water Temperature and Dissolved Organic Matter on Ozone and OH Radical Reaction Kinetics During Ozonation of a Lake Water. *Ozone: Science & Engineering* 2016, 38, 100-114. DOI: 10.1080/01919512.2015.1079120.
37. Yury Gorbanev, Deborah O'Connell, Victor Chechik. Non-Thermal Plasma in Contact with Water: The Origin of Species. *Chemistry - A European Journal* 2016, 22 (10) , 3496-3505. DOI: 10.1002/chem.201503771.
38. Arnold N. Onyango. Alternatives to the 'water oxidation pathway' of biological ozone formation. *Journal of Chemical Biology* 2016, 9, 1-8. DOI: 10.1007/s12154-015-0140-6.
39. Fei Qi, Wei Chu, Bingbing. u. Comparison of phenacetin degradation in aqueous solutions by catalytic ozonation with CuFe₂O₄ and its precursor: Surface properties, intermediates and reaction mechanisms. *Chemical Engineering Journal* 2016, 284, 28-36. DOI: 10.1016/j.cej.2015.07.095.
40. Erika Reisz, Sergej Naumov, Winfried Schmidt, Clemens von Sonntag. Reaction of Ozone with Ag(I)—Mechanistic Considerations. *Ozone: Science & Engineering* 2015, 37, 393-404. DOI: 10.1080/01919512.2015.1041583.
41. Huijiao Wang, Shi Yuan, Juhong Zhan, Yujue Wang, Gang Yu, Shubo Deng, Jun Huang, Bin Wang. Mechanisms of enhanced total organic carbon elimination from oxalic acid solutions by electro-peroxone process. *Water Research* 2015, 80, 20-29. DOI: 10.1016/j.watres.2015.05.024.
42. Jesús Ferre-Aracil, Salvador C. Cardona, Javier Navarro-Laboulais. Kinetic Study of Ozone Decay in Homogeneous Phosphate-Buffered Medium. *Ozone: Science & Engineering* 2015, 37, 330-342. DOI: 10.1080/01919512.2014.998756.
43. Fei Qi, Bingbing. u, Wei Chu. Heterogeneous catalytic ozonation of phenacetin in water using magnetic spinel ferrite as catalyst: Comparison of surface property and efficiency. *Journal of Molecular Catalysis A: Chemical* 2015, 396, 164-173. DOI: 10.1016/j.molcata.2014.10.001.
44. Erika Reisz, Alexandra Fischbacher, Sergej Naumov, Clemens von Sonntag, Torsten C. Schmidt. Hydride Transfer: A Dominating Reaction of Ozone with Tertiary Butanol and Formate Ion in Aqueous Solution. *Ozone: Science & Engineering* 2014, 36, 532-539. DOI: 10.1080/01919512.2014.891436.
45. Sergey L. Khursan. Peroxide intermediates of oxidation processes: Organic trioxides. 2014,,, 1-72. DOI: 10.1002/9780470682531.pat0874.

46. S. Reed Plimpton, Mark Gołkowski, Deborah G. Mitchell, Chad Austin, Sandra S. Eaton, Gareth R. Eaton, Czesław Gołkowski, Martin Voskuil. Remote delivery of hydroxyl radicals via secondary chemistry of a nonthermal plasma effluent. *Biotechnology and Bioengineering* 2013, 110 (10.1002/bit.v110.7) , 1936-1944. DOI: 10.1002/bit.24853.
47. Zdenko Machala, Barbora Tarabova, Karol Hensel, Eva Spetlikova, Libusa Sikurova, Petr Lukes. Formation of ROS and RNS in Water Electro-Sprayed through Transient Spark Discharge in Air and their Bactericidal Effects. *Plasma Processes and Polymers* 2013, 10 (10.1002/ppap.v10.7) , 649-659. DOI: 10.1002/ppap.201200113.
48. Josep M. Anglada, Miquel Torrent-Sucarrat, Manuel F. Ruiz-Lopez, Marilia Martins-Costa. Is the HO₄⁻ Anion a Key Species in the Aqueous-Phase Decomposition of Ozone?. *Chemistry - A European Journal* 2012, 18 (10.1002/chem.v18.42) , 13435-13445. DOI: 10.1002/chem.201201991.
49. Marek Golkowski, Czesław Golkowski, Jori Leszczynski, S. Reed Plimpton, Piotr Maslowski, Aleksandra Foltynowicz, Jun Ye, Bruce McCollister. Hydrogen-Peroxide-Enhanced Nonthermal Plasma Effluent for Biomedical Applications. *IEEE Transactions on Plasma Science* 2012, 40, 1984-1991. DOI: 10.1109/TPS.2012.2200910.
50. D. Gardoni, A. Vailati, R. Canziani. Decay of Ozone in Water: A Review. *Ozone: Science & Engineering* 2012, 34, 233-242. DOI: 10.1080/01919512.2012.686354.
51. Sergej Naumov, Clemens von Sonntag. The reaction of •OH with O₂, the decay of O₃^{•-} and the pK_a of HO₃[•] - interrelated questions in aqueous free-radical chemistry. *Journal of Physical Organic Chemistry* 2011, 24 (10.1002/poc.v24.7) , 600-602. DOI: 10.1002/poc.1812.
52. Rossano Amadelli, Luca Samiolo, Achille De Battisti, Alexander B. Velichenko. Electro-oxidation of Some Phenolic Compounds by Electrogenerated O₃ and by Direct Electrolysis at PbO₂ Anodes. *Journal of The Electrochemical Society* 2011, 158, P87. DOI: 10.1149/1.3589913.
53. Pentax, White Papers
54. ASP, White Papers
55. Olympus, White Papers